



DWYER PRECISION® PRODUCTS INC.

266 N. 20th Street
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

| BREATHCALL® | | Retail Price List | | Sept 1, 2011 |
|--------------------|---------------------------|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No. | Item Photo. | Price. | Systems Used In. | Electrical Characteristics |
| 11001-09' | New Multi Singles | 136.15 | IN-LINE Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11201-09' | New Multi Singles | 136.15 | IN-LINE Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11050-09' | Classic Model Singles | 146.00 | Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in... | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11250-09' | Classic Model Singles | 146.00 | Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in... | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11800-09' | Classic Model Singles | 137.25 | Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800 | Normally closed contact, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11900-09' | Classic Model Singles | 137.25 | Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800 | Normally closed contact, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★







NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

DWYER PRECISION® PRODUCTS INC.

266 N. 20th Street
Jacksonville Beach, Florida 32250







All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

| BREATHCALL® | | Retail Price List | | Sept 1, 2011 |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No. | Item Photo. | Price. | Systems Used In. | Electrical Characteristics |
| 11001-09'D | New Multi Duals  Classic Model Nurse TV | 182.45 | Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11201-09'D | New Multi Duals  Classic Model Nurse TV | 182.45 | Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11050-09'D | Classic Model Duals  Nurse TV | 200.60 | Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in... | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11250-09'D | Classic Model Duals  Nurse TV | 200.60 | Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in... | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11800-09'D | Classic Model Duals  Nurse TV | 190.80 | Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800 | Normally closed contact, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11900-09'D | Classic Model Duals  Nurse TV | 190.80 | Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800 | Normally closed contact, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

| BREATHCALL® | | Retail Price List | | Sept 1, 2011 |
|--------------------|--------------------------------------------------------------------------------------------------------------|--------------------------|------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No. | Item Photo. | Price. | Systems Used In. | Electrical Characteristics |
| 11001-EK8-09' | Classic Model Singles  | 153.50 | Ektacom - 8 - Pin Octual In-line with Keyway | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-EK15-09' | Classic Model Singles  | 155.60 | Fisher Berkely - 15 - Pin Amphenol NMC / Lifeline In-Line Model | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11701-AZ6-09' | Classic Model Singles  | 171.90 | American Zettler - 6 - Contact (Green Spool) In-Line Model | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-DN8-09'JR | Classic Model Singles  | 156.40 | 8 Pin Din Plug Jeron System Standard Nurse Call with resistor | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-AZ9-09' | Classic Model Singles  | 159.20 | American Zettler - 9 - Contact In-Line Model | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11800-AZ9-09' | | 189.10 | Bunting Style NC (closed) Must supply pin view | |
| 11001-SY20-09' | Classic Model Singles  | 153.50 | HCE- Sylvania / Phillips - 20 - Pin Plug In-Line Model | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |







★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

BREATHCALL®

Retail Price List

Sept 1, 2011

| Catalog No. | Item Photo. | Price. | Systems Used In. | Electrical Characteristics |
|----------------|--------------------------------------------------------------------------------------------------------------|--------|--------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 11001-AK7-09' | Classic Model Singles  | 150.70 | 7 - Pin Din Plug Ackerman System Standard Nurse Call | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-CD7-09' | Classic Model Singles  | 150.70 | 7 - Pin Din Plug Condigi Model Standard Nurse Call | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-OS5-09' | Classic Model Singles  | 150.70 | 5 - Pin Din Plug Ohmeda Schmelter Standard Nurse Call | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-DN8-09'R | Classic Model Singles  | 150.70 | 8 - Pin Din Plug Rauland Borg System Standard Nurse Call | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-DN8-09'J | Classic Model Singles  | 150.70 | 8 - Pin Din Plug Jerom Systems Standard Nurse Call | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11800-AK7-09' |  | 183.80 | 7 - Pin Din Plug Ackermann System (German) Standard Nurse Call Bunting Style / NC | Normally closed contact, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |







★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

BREATHCALL®

Retail Price List

Sept 1, 2011

| Catalog No. | Item Photo. | Price. | Systems Used In. | Electrical Characteristics |
|---------------|--------------------------------------------------------------------------------------------------------------|--------|--------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 11050-AU5-09' | Classic Model Singles  | 165.95 | Auth - 5 - Pin Plug Maintain Contact | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11050-FA5-09' | Classic Model Singles  | 165.95 | Faraday - 5 - Pin Plug Maintain Contact | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11050-ED5-09' | Classic Model Singles  | 165.95 | Edwards - 5 - Pin Plug Maintain Contact | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11050-C05-09' | Classic Model Singles  | 165.95 | Couch - 5 - Pin Plug Maintain Contact | Normally open contacts, blowing into soft straw mouthpiece makes contact that is maintained until reset button (on top of plug at wall station) is pushed down (depressed). No system relay is required. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-ED3-09' | Classic Model Singles  | 153.50 | Edwards - 3 - Prong Plug Maintain Contact | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |
| 11001-WE8-09' | Classic Model Singles  | 153.50 | Intego / Wescom / Tektone / Simplex 8 pin connector | Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC |

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)