

DWYER PRECISION® PRODUCTS INC.

266 N. 20th Street
Jacksonville Beach, Florida 32250







All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

HEADCALL®

Retail Price List

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
15007-06'	 Singles	40.90	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing pad makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
15007-08'		46.40		
15007-10'		51.90		
15007-12'		57.40		
15207-12'	 Singles	41.50	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing pad makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
15207-08'		47.00		
15207-10'		52.50		
15207-12'		58.00		
15201-06'	 Singles	38.00	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing pad makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
15201-08'		43.50		
15201-10'		49.00		
15201-12'		54.50		
15027-06'	 Duals	44.20	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing pad makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
15027-08'		49.70		
15027-10'		55.20		
15027-12'		60.70		
15225-06'	 Duals	43.50	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing pad makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
15225-08'		50.70		
15225-10'		57.90		
15225-12'		65.00		
15227-06'	 Duals	45.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing pad makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
15227-08'		52.60		
15227-10'		59.80		
15227-12'		67.00		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)