













<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13001A-06'		47.15	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13001A-08'		49.35		
13001A-10'		52.65		
13001A-12'		53.80		
13007A-06'		50.45	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007A-08'		52.65		
13007A-10'		53.80		
13007A-12'		57.05		
13025A-06'		72.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13025A-08'		76.85		
13025A-10'		81.20		
13027A-06'		85.65	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027A-08'		90.20		
13027A-10'		94.55		
13050-06'		61.15	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13050-08'		63.55		
13050-10'		65.25		
13050-12'		68.30		
13075-06'		84.75	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13075-08'		89.55		
13075-10'		94.40		
13075-12'		99.05		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★







NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13900-06'	 <p>Classic Model Singles</p>	51.45	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
13900-08'		53.15		
13900-10'		55.35		
13900-12'		57.55		
13925-06'	 <p>Classic Model Duals</p>	75.30	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
13925-08'		79.65		
13925-10'		84.15		
13925-12'		88.20		
13201B-06'	 <p>Singles</p>	46.75	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13201B-08'		48.90		
13201B-10'		52.15		
13201B-12'		54.35		
13207B-06'	 <p>Singles</p>	50.00	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13207B-08'		51.25		
13207B-10'		53.55		
13207B-12'		55.40		
13225B-06'	 <p>Duals</p>	70.30	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13225B-08'		74.60		
13225B-10'		78.85		
13250-06'	 <p>Classic Model Singles</p>	59.35	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13250-08'		61.70		
13250-10'		64.00		
13250-12'		66.35		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)



<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13050-AU5-06'	 <p>Classic Model Singles</p>	78.60	Auth - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13050-AU5-08'		80.90		
13050-AU5-10'		83.20		
13050-AU5-12'		85.45		
13075-AU5-06'	 <p>Classic Model Duals</p>	91.50	Auth - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13075-AU5-08'		95.50		
13075-AU5-10'		100.05		
13075-AU5-12'		104.30		
13050-ED5-06'	 <p>Classic Model Singles</p>	78.60	Edwards - 5 - Prong Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13050-ED5-08'		80.90		
13050-ED5-10'		83.20		
13050-ED5-12'		85.45		
13075-ED5-06'	 <p>Classic Model Duals</p>	91.50	Edwards - 5 - Prong Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13075-ED5-08'		95.80		
13075-ED5-10'		100.05		
13075-ED5-12'		104.30		
13050-FA5-06'	 <p>Classic Model Singles</p>	78.60	Faraday - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13050-FA5-08'		80.90		
13050-FA5-10'		83.20		
13050-FA5-12'		85.45		
13075-FA5-06'	 <p>Classic Model Duals</p>	91.50	Faraday - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13075-FA5-08'		95.80		
13075-FA5-10'		100.05		
13075-FA5-12'		104.30		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13050-CO5-06'	 <p>Classic Model Singles</p>	78.60	Couch - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13050-CO5-08'		80.90		
13050-CO5-10'		83.20		
13050-CO5-12'		85.45		
13075-CO5-06'	 <p>Classic Model Duals</p>	91.50	Couch - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13075-CO5-08'		95.50		
13075-CO5-10'		100.05		
13075-CO5-12'		104.30		
13007-ED3-06'	 <p>Classic Model Singles</p>	67.55	Edwards - 3 - Prong Plug Inline (Lifeline model)	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-ED3-08'		69.85		
13007-ED3-10'		72.10		
13007-ED3-12'		74.70		
13027-ED3-06'	 <p>Classic Model Duals</p>	91.50	Edwards - 3 - Prong Plug Inline (Lifeline model)	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-ED3-08'		95.50		
13027-ED3-10'		100.05		
13027-ED3-12'		104.30		
13007-WE8-06'	 <p>Classic Model Singles</p>	66.40	Intego / Wescom / Tektone / Simplex 8 pin connector	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-WE8-08'		68.65		
13007-WE8-10'		70.95		
13007-WE8-12'		73.25		
13027-WE8-06'	 <p>Classic Model Duals</p>	87.80	Intego / Wescom / Tektone / Simplex 8 pin connector	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-WE8-08'		92.50		
13027-WE8-10'		95.25		
13027-WE8-12'		98.75		







★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

**PADCALL®**

**Retail Price List**

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13007-EK8-06'	 <p>Classic Model Singles</p>	68.65	Ektacom - 8 - Pin Octual In-Line Model with Keyway	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-EK8-08'		70.95		
13007-EK8-10'		73.25		
13007-EK8-12'		76.55		
13027-EK8-06'	 <p>Classic Model Duals</p>	89.70	Ektacom - 8 - Pin Octual In-Line Model with Keyway	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-EK8-08'		94.50		
13027-EK8-10'		99.20		
13027-EK8-12'		103.95		
13007-EK15-06'	 <p>Classic Model Singles</p>	80.00	Fisher Berkely - 15 - Pin Amphenol NMC / Lifeline In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-EK15-08'		82.15		
13007-EK15-10'		84.40		
13007-EK15-12'		86.60		
13027-EK15-06'	 <p>Classic Model Duals</p>	91.50	Fisher Berkely - 15 - Pin Amphenol NMC / Lifeline In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-EK15-08'		95.50		
13027-EK15-10'		100.05		
13027-EK15-12'		104.30		
13707-AZ6-06'	 <p>Classic Model Singles</p>	91.50	American Zettler - 6 - Contact (Green Spool) In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13707-AZ6-08'		95.50		
13707-AZ6-10'		100.05		
13707-AZ6-12'		104.30		
13707-AZ6-06'	 <p>Classic Model Singles</p>		<u>DISCONTINUED</u>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13707-AZ6-06'			American Zettler - 4 - Contact (White Spool)	
13707-AZ6-10'			In-Line Model	
13707-AZ6-12'				

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250







All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13727-AZ6-06'	 Classic Model Duals	109.95	American Zettler - 6 - Contact (Green Spool) In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13727-AZ6-08'		114.80		
13727-AZ6-10'		119.70		
13727-AZ6-12'		124.75		
13727-AZ6-06'	 Classic Model Duals		<u>DISCONTINUED</u>  American Zettler - 4 - Contact (White Spool) In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13727-AZ6-08'				
13727-AZ6-10'				
13727-AZ6-12'				
13007-AZ9-06'	 Classic Model Singles	73.85	American Zettler - 9- Contact In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-AZ9-08'		76.20		
13007-AZ9-10'		78.45		
13007-AZ9-12'		80.80		
13027-AZ9-06'	 Classic Model Duals	91.50	American Zettler - 9- Contact In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-AZ9-08'		95.50		
13027-AZ9-10'		100.05		
13027-AZ9-12'		104.30		
13007-SY20-06'	 Classic Model Singles	66.40	HCE- Sylvania / Phillips - 20 - Pin Plug In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-SY20-08'		68.65		
13007-SY20-10'		70.95		
13007-SY20-12'		73.25		
13027-SY20-06'	 Classic Model Duals	89.70	HCE- Sylvania / Phillips - 20 - Pin Plug In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-SY20-08'		94.50		
13027-SY20-10'		99.20		
13027-SY20-12'		104.00		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13050-SY20-06'	 <p>Classic Model Singles</p>	76.20	HCE- Sylvania / Phillips - 20 - Pin Maintain Contract (Right Angle Only)	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13050-SY20-08'		78.35		
13050-SY20-10'		80.60		
13050-SY20-12'		82.85		
13075-SY20-06'	 <p>Classic Model Duals</p>	92.25	HCE- Sylvania / Phillips - 20 - Pin Maintain Contract (Right Angle Only)	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
13075-SY20-08'		97.05		
13075-SY20-10'		101.80		
13075-SY20-12'		106.50		
13001-AU-LT-06'	 <p>Classic Model Singles</p>		<i>DISCONTINUED</i>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13001-AU-LT-08'				
13001-AU-LT-10'				
13001-AU-LT-12'				
13001-AU-RT-06'	 <p>Classic Model Singles</p>		<i>DISCONTINUED</i>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13001-AU-RT-08'				
13001-AU-RT-10'				
13001-AU-RT-12'				
13007-DN8-06'R	 <p>Classic Model Singles</p>	63.85	8 - Pin Din Plug Rauland Systems Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-DN8-08'R		66.10		
13007-DN8-10'R		68.30		
13007-DN8-12'R		70.55		
13007-DN8-06'J	 <p>Classic Model Singles</p>	63.85	8 - Pin Din Plug Jeron Systems Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-DN8-08'J		66.10		
13007-DN8-10'J		68.30		
13007-DN8-12'J		70.55		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

*All Nurse Call cords are suitable for use in oxygen atmosphere.*

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13007-DN8-06'JR	 <p>Classic Model Singles</p>	66.65	8 - Pin Din Plug Jeron Systems Resistor Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-DN8-08'JR		69.95		
13007-DN8-10'JR		72.30		
13007-DN8-12'JR		74.50		
13027-DN8-06'-R	 <p>Classic Model Duals</p>	89.45	8 - Pin Din Plug Rauland Borg System Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-DN8-08'-R		94.25		
13027-DN8-10'-R		98.95		
13027-DN8-12'-R		103.70		
13027-DN8-06'-J	 <p>Classic Model Duals</p>	89.45	8 - Pin Din Plug Jeron Systems Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-DN8-08'-J		94.25		
13027-DN8-10'-J		98.95		
13027-DN8-12'-J		103.70		
13027-DN8-06'JR	 <p>Classic Model Duals</p>	92.60	8 - Pin Din Plug Jeron Systems Resistor Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-DN8-08'JR		100.65		
13027-DN8-10'JR		103.20		
13027-DN8-12'JR		106.85		
13800-DN8-06'	 <p>Classic Model Singles</p>	69.35	8 - Pin Din Plug Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Right angle only..	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
13800-DN8-08'		71.60		
13800-DN8-10'		73.85		
13800-DN8-12'		76.10		
13827-DN8-06'	 <p>Classic Model Duals</p>	93.90	8 - Pin Din Plug Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Right angle only..	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
13827-DN8-08'		98.35		
13827-DN8-10'		100.60		
13827-DN8-12'		102.80		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PADCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
13007-AK7-06'	 Classic Model Singles	64.10	7 - Pin Din Plug Ackerman System Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-AK7-08'		66.30		
13007-AK7-10'		68.45		
13007-AK7-12'		70.65		
13027-AK7-06'	 Classic Model Duals	90.25	7 - Pin Din Plug Ackerman System Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-AK7-08'		94.55		
13027-AK7-10'		98.90		
13027-AK7-12'		104.35		
13007-CD7-06'	 Classic Model Singles	64.10	7 - Pin Din Plug Condigi Model Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-CD7-08'		66.30		
13007-CD7-10'		68.45		
13007-CD7-12'		70.65		
13027-CD7-06'	 Classic Model Duals	90.25	7 - Pin Din Plug Condigi Model Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-CD7-08'		94.55		
13027-CD7-10'		98.90		
13027-CD7-12'		104.35		
13007-OS5-06'	 Classic Model Singles	64.10	5 - Pin Din Ohmeda Schmelter Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13007-OS5-08'		66.30		
13007-OS5-10'		68.45		
13007-OS5-12'		70.65		
13027-OS5-06'	 Classic Model Duals	90.25	5 - Pin Din Ohmeda Schmelter Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
13027-OS5-08'		94.55		
13027-OS5-10'		99.00		
13027-OS5-12'		104.35		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)