









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>	
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics	
10001A-06'	 Singles	32.30	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC	
10001A-08'		34.45			
10001A-10'		37.50			
10007A-06'	 Singles	38.35	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC	
10007A-08'		38.90			
10007A-10'		39.55			
10025A-06'	 Duals	46.45	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC	
10025A-08'		51.20			
10027A-06'	 Duals	49.85	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC	
10027A-08'		54.95			
10050-06'	 Classic Model Singles	46.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC	
10050-08'		48.80			
10050-10'		51.45			
10050-12'		53.90			
10075-06'	 Classic Model Duals	60.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC	
10075-08'		65.20			
10075-10'		69.45			
10075-12'		75.45			

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)










**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

		<h2 style="text-align: center;">Retail Price List</h2>		<h2 style="margin: 0;">Sept 1, 2011</h2>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10900-06'	 <p>Classic Model Singles</p>	36.20	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10900-08'		41.25		
10900-10'		44.30		
10900-12'		46.75		
10925-06'	 <p>Classic Model Duals</p>	62.10	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10925-08'		67.00		
10925-10'		71.80		
10925-12'		79.00		
10201B-06'	 <p>Singles</p>	32.35	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10201B-08'		34.50		
10201B-10'		37.55		
10207B-06'	 <p>Singles</p>	38.95	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 12' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10207B-08'		39.65		
10207B-10'		40.20		
10225B-06'	 <p>Duals</p>	46.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10225B-08'		51.20		
10227B-06'	 <p>Duals</p>	60.60	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.  <b>Note: Unit not available over 10' foot</b>	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10227B-08'		65.40		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10250-06'	 Classic Model Singles	46.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10250-08'		48.80		
10250-10'		50.55		
10250-12'		53.90		
10275-06'	 Classic Model Duals	60.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10275-08'		65.20		
10275-10'		69.80		
10275-12'		75.45		
10800-06'	 Classic Model Singles	38.40	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10800-08'		41.25		
10800-10'		44.30		
10800-12'		46.75		
10825-06'	 Classic Model Duals	62.10	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10825-08'		67.00		
10825-10'		71.80		
10825-12'		77.50		
10007-06'D	 Classic Model with TV Duals  Nurse TV	93.95	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. with TV	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-08'D		95.80		
10007-10'D		97.75		
10007-12'D		99.75		
10207-06'D	 Classic Model with TV Duals  Nurse TV	93.95	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std. with TV	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10207-08'D		95.80		
10207-10'D		97.75		
10207-12'D		99.75		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10007-DN8-06'D	 <p>Duals</p> <p>Nurse TV</p>	108.65	8-pin Din Plug Rauland and Jeron Systems Standard Nurse Call. with TV	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-DN8-08'D		110.00		
10007-DN8-10'D		112.25		
10007-DN8-12'D		121.00		
10007-WE8-06'D	 <p>Duals</p> <p>Nurse TV</p>	108.65	Intego / Wescom / Tektone / Simplex 8 pin connector with TV	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-WE8-08'D		110.00		
10007-WE8-10'D		112.25		
10007-WE8-12'D		121.00		
10050-AU5-06'	 <p>Classic Model</p> <p>Singles</p>	76.80	Auth - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10050-AU5-08'		79.15		
10050-AU5-10'		81.40		
10050-AU5-12'		83.60		
10075-AU5-06'	 <p>Classic Model</p> <p>Duals</p>	94.45	Auth - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10075-AU5-08'		97.70		
10075-AU5-10'		101.10		
10075-AU5-12'		103.45		
10050-ED5-06'	 <p>Classic Model</p> <p>Singles</p>	76.80	Edwards - 5 - Prong Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10050-ED5-08'		79.15		
10050-ED5-10'		81.40		
10050-ED5-12'		83.60		
10075-ED5-06'	 <p>Classic Model</p> <p>Duals</p>	94.45	Edwards - 5 - Prong Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10075-ED5--08'		97.70		
10075-ED5-10'		101.10		
10075-ED5-12'		103.45		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)



**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10050-C05-06'	<p>Singles</p>	76.80	Couch - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10050-C05-08'		79.15		
10050-C05-10'		81.40		
10050-C05-12'		83.60		
10075-C05-06'	<p>Duals</p>	84.00	Couch - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10075-C05-08'		88.10		
10075-C05-10'		92.70		
10075-C05-12'		97.20		
10007-ED3-06'	<p>Singles</p>	58.15	Edwards - 3 - Prong Plug In-Line (Lifeline model)	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-ED3-08'		59.85		
10007-ED3-10'		62.17		
10007-ED3-12'		64.40		
10027-ED3-06'	<p>Duals</p>	80.20	Edwards - 3 - Prong Plug In-Line (Lifeline model)	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-ED3-08'		85.90		
10027-ED3-10'		90.40		
10027-ED3-12'		95.80		
10007-WE8-06'	<p>Singles</p>	57.05	Intego / Wescom / Tektone / Simplex 8 pin connector	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-WE8-08'		58.75		
10007-WE8-10'		62.15		
10007-WE8-12'		64.40		
10027-WE8-06'	<p>Duals</p>	66.85	Intego / Wescom / Tektone / Simplex 8 pin connector	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-WE8-08'		70.55		
10027-WE8-10'		75.60		
10027-WE8-12'		82.95		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10800-WE8-06'	<p>Singles</p>	65.70	Intego / Wescom / Tektone / Simplex 8 pin connector Bunting Style NC/ (right angle only).	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10800-WE8-08'		68.00		
10800-WE8-10'		70.90		
10800-WE8-12'		73.15		
10827-WE8-06'	<p>Duals</p>	76.15	Intego / Wescom / Tektone / Simplex 8 pin connector Bunting Style NC/ (right angle only).	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10827-WE8-08'		81.15		
10827-WE8-10'		86.75		
10827-WE8-12'		91.80		
10007-EK8-06'	<p>Singles</p>	56.65	Ektacom - 8 - Pin Octual In-Line with Keyway	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-EK8-08'		58.95		
10007-EK8-10'		61.30		
10007-EK8-12'		63.60		
10027-EK8-06'	<p>Duals</p>	81.05	Ektacom - 8 - Pin Octual In-Line with Keyway	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-EK8-08'		85.90		
10027-EK8-10'		90.70		
10027-EK8-12'		96.10		
10007-EK15-06'	<p>Singles</p>	70.05	Fisher Berkely - 15 - Pin Amphenol MMC / Lifeline In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-EK15-08'		72.30		
10007-EK15-10'		74.55		
10007-EK15-12'		76.80		
10027-EK15-06'	<p>Duals</p>	81.90	Fisher Berkely - 15 - Pin Amphenol MMC / Lifeline In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-EK15-08'		86.80		
10027-EK15-10'		91.55		
10027-EK15-12'		97.50		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10707-AZ6-06'	 <p>Classic Model Singles</p>	83.30	American Zettler - 6 - Contact (Green Spool) In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10707-AZ6-08'		85.80		
10707-AZ6-10'		88.20		
10707-AZ6-12'		90.60		
10727-AZ6-06'	 <p>Classic Model Duals</p>	88.40	American Zettler - 6 - Contact (Green Spool) In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10727-AZ6-08'		93.70		
10727-AZ6-10'		99.00		
10727-AZ6-12'		105.40		
10007-AZ9-06'	 <p>Classic Model Singles</p>	66.75	American Zettler - 9 - Contact In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-AZ9-08'		69.10		
10007-AZ9-10'		71.50		
10007-AZ9-12'		73.95		
10027-AZ9-06'	 <p>Classic Model Duals</p>	75.70	American Zettler - 9 - Contact In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-AZ9-08'		80.85		
10027-AZ9-10'		86.05		
10027-AZ9-12'		92.15		
10001-06'	 <p>Classic Model Singles</p>	33.85	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10001-08'		36.00		
10001-10'		39.05		
10001-12'		41.45		
10201-06'	 <p>Classic Model Singles</p>	33.90	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10201-08'		36.05		
10201-10'		39.10		
10201-12'		41.50		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit). Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10800-AZ9-06'	 Classic Model Singles	72.50	American Zettler - 9 - Contact Amphenol (Right Angle Only) Bunting Style NC/	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10800-AZ9-08'		75.00		
10800-AZ9-10'		77.50		
10800-AZ9-12'		79.75		
10827-AZ9-06'	 Classic Model Duals	86.35	American Zettler - 9 - Contact Amphenol (Right Angle Only) Bunting Style NC/	Normally closed contact, pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed.  Rated 1/3 amps 24 volts AC/DC. Suitable for use in oxygen atmosphere.
10827-AZ9-08'		88.68		
10827-AZ9-10'		91.05		
10827-AZ9-12'		93.30		
10007-SY20-06'	 Classic Model Singles	56.45	HCE- Sylvania / Phillips - 20 - Pin Plug In-Line	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-SY20-08'		58.70		
10007-SY20-10'		61.00		
10007-SY20-12'		63.20		
10027-SY20-06'	 Classic Model Duals	72.85	HCE- Sylvania / Phillips - 20 - Pin Plug In-Line	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-SY20-08'		77.35		
10027-SY20-10'		81.75		
10027-SY20-12'		87.00		
10050-SY20-06'	 Classic Model Singles	70.80	HCE- Sylvania / Phillips - 20 - Pin Plug Maintain Contact (Right Angle Only)	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10050-SY20-08'		72.00		
10050-SY20-10'		74.30		
10050-SY20-12'		79.10		
10075-SY20-06'	 Classic Model Duals	85.75	HCE- Sylvania / Phillips - 20 - Pin Plug Maintain Contact (Right Angle Only)	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10075-SY20-08'		90.25		
10075-SY20-10'		93.90		
10075-SY20-12'		98.55		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10001-AU-LT-06'	 <p>Classic Model Singles</p>		<p><i>DISCONTINUED</i></p> <p>Standard Auth 1/4" Phone Plug Left Side 5/8" center</p>	<p>Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC</p>
10001-AU-LT-08'				
10001-AU-LT-10'				
10001-AU-LT-12'				
10001-AU-RT-06'	 <p>Classic Model Singles</p>		<p><i>DISCONTINUED</i></p> <p>Standard Auth 1/4" Phone Plug Left Side 5/8" center</p>	<p>Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC</p>
10001-AU-RT-08'				
10001-AU-RT-10'				
10001-AU-RT-12'				
10007-AK7-06'	 <p>Singles</p>	54.60	<p>7 - Pin Din Plug Ackerman System Standard Nurse Call</p>	<p>Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC</p>
10007-AK7-08'		56.70		
10007-AK7-10'		58.80		
10007-AK7-12'		61.65		
10027-AK7-06'	 <p>Duals</p>	67.05	<p>7 - Pin Din Plug Ackerman System Standard Nurse Call</p>	<p>Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC</p>
10027-AK7-08'		72.10		
10027-AK7-10'		77.00		
10027-AK7-12'		82.95		
10007-DN8-06'J	 <p>Singles</p>	54.60	<p>8 - Pin Din Plug Jeron Systems Standard Nurse Call</p>	<p>Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC</p>
10007-DN8-08'J		56.70		
10007-DN8-10'J		58.80		
10007-DN8-12'J		61.65		
10027-DN8-06'J	 <p>Duals</p>	67.05	<p>8 - Pin Din Plug Jeron Systems Standard Nurse Call</p>	<p>Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC</p>
10027-DN8-08'J		72.10		
10027-DN8-10'J		77.00		
10027-DN8-12'J		82.95		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









**DWYER PRECISION® PRODUCTS INC.**

266 N. 20th Street  
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

<b>PRESSCALL®</b>		<b>Retail Price List</b>		<b>Sept 1, 2011</b>
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
10007-DN8-06'-R	 Classic Model Singles	54.60	8 - Pin Din Plug Rauland Borg System Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-DN8-08'-R		56.70		
10007-DN8-10'-R		58.80		
10007-DN8-12'-R		61.65		
10027-DN8-06'-R	 Classic Model Duals	67.05	8 - Pin Din Plug Rauland Borg System Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-DN8-08'-R		72.10		
10027-DN8-10'-R		77.00		
10027-DN8-12'-R		82.95		
10007-DN8-06'-JR	 Classic Model Singles	60.50	8 - Pin Din Plug Jeron System Resistor Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10007-DN8-08'-JR		62.70		
10007-DN8-10'-JR		65.30		
10007-DN8-12'-JR		67.50		
10027-DN8-06'-JR	 Classic Model Duals	72.85	8 - Pin Din Plug Jeron System Resistor Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
10027-DN8-08'-JR		77.90		
10027-DN8-10'-JR		82.85		
10027-DN8-12'-JR		88.70		
10050-FA5-06'	 Classic Model Singles	77.55	Faraday - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10050-FA5-08'		79.90		
10050-FA5-10'		82.20		
10050-FA5-12'		84.45		
10075-FA5-06'	 Classic Model Duals	88.90	Faraday - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC  Suitable for use in oxygen atmosphere.
10075-FA5-08'		93.25		
10075-FA5-10'		95.15		
10075-FA5-12'		98.95		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).  
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)