

REPLACEMENT PLUG KITS

Retail Price List

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
113-10207	Classic Model PK 	32.60	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
113-10250	Classic Model PK 	30.60	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-10201	Classic Model PK 	29.55	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
113-10800	Classic Model PK 	33.65	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact - Pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-10900	Classic Model PK 	33.65	Bunting Sterl System for lighted cord stations. (Presscall bulb is NOT lighted). Electrical plug has plastic tip, brass ring and sleeve model #10900. (ring and sleeve contact) Tip and sleeve contact Model #10800	Normally closed contact - Pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-DN8	Classic Model PK 	33.55	(MUST STATE WHICH SYSTEM) 8 - Pin Din Plug Rauland / Jeron Systems Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC







★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

REPLACEMENT PLUG KITS

Retail Price List

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
113-AU5	Classic Model PK 	34.65	Auth - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-FA5	Classic Model PK 	35.70	Faraday - 5 - Pin Plug Mainain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-ED5	Classic Model PK 	35.70	Edwards - 5 - Prong Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-C05	Classic Model PK 	35.70	Couch - 5 - Prong Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
113-ED3	Classic Model PK 	32.60	Edwards - 3 - Prong Plug In-Line (Lifeline model)	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC
113-WE8	Classic Model Singles 	29.55	Intego / Wescom / Tektone / Simplex 8 pin connector	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC







★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

REPLACEMENT PLUG KITS

Retail Price List

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
113-EK8	 <p>Classic Model PK</p>	33.65	Ektacom - 8 - Pin Octual In-Line Model with Keyway	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC
113-EK15	 <p>Classic Model PK</p>	36.70	Fisher Berkely - 15 - Pin Amphenol NMC / Lifeline In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC
113-AZ9	 <p>Classic Model PK</p>	36.70	American Zettler - 9 - Contact Amphenol with Pigtail In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC
113-SY20	 <p>Classic Model PK</p>	35.70	HCE- Sylvania / Phillips - 20 - Pin Plug In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC
113-AZ6	 <p>Classic Model PK</p>	48.35	American Zettler - 6 - Contact (Green Spool) In-Line Model	Momentary contact type. Normally open contacts, blowing into soft straw mouthpiece makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug is rated 1/3 amps 24 volt AC/DC
113-DN8B	 <p>Classic Model PK</p>	35.35	Bunting 8 Pin DIN for Jeron & Rauland Systems (Right Angle Only)	Normally closed contact - Pressing bulb breaks contact momentarily and system switches over to call circuit and holds until reset button at bed station or central nurse station is depressed. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.






★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)

**REPLACEMENT FILTER KITS
and MULTI PURPOSE CLAMP**

Retail Price List

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics		
113-21001	 <p>1 Unit</p>	32.60	<p><u>CLEAR</u></p> <p>Accessory Packages are used with All Breathcall Dual Breathcall (Nurse/TV) and all Mini Breathcall Products.</p>	<p>All Accessory Packages come with 12 Filter Assemblys with softstraws & Porous Plastic Filters. All Filter assemblys should be replaced after third cleaning.</p>		
	Filter Kit					
113-21002	 <p>1 Unit</p>					
	Filter Kit					
113-21003	 <p>1 Unit</p>		<p><u>CLEAR</u></p> <p>Accessory Packages are used with All Breathcall Dual Breathcall (Nurse/TV) and all Mini Breathcall Products.</p>	<p>All Accessory Packages come with 12 Filter Assemblys with softstraws & Porous Plastic Filters. All Filter assemblys should be replaced after third cleaning.</p>		
	Filter Kit					
113-11020	 <p>1 Unit</p>	22.45			<p>Used on all Breathcall and Dual Breathcall Units, also can be used for lightweight clamping for all sorts of Hobbies where an extra hand could be used.</p>	<p>Multi-Purpose Assembly mounting Clamp, 1 unit "C" Shaped for great gripping strength (med:) Dimension Approx 3-5 1/8" x 2" at the bottom and 1-7/8" at the top. Cutout area 1-1/2" by 2-1/4".</p>
	Clamp					
113-21004		2.50	<p>CLEAR Flex Tubes</p> <p>Used on Dual "Y" Nurse/TV as adjustable extension for filter assembly, slip fit over "Y" ends and insert filter at top per photo.</p> <p>(Suggested use for ridged "Y" connections only)</p>	<p>Expandable - flexible filter tube extension adapter.</p> <p>DISPOSABLE</p> <p>Expands from 2-3/4" to 5-3/8" in length and fully adjustable to suit your patient comfort and needs. 12-Units to a pack (filters not included).</p>		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★