

DWYER PRECISION® PRODUCTS INC.

266 N. 20th Street
Jacksonville Beach, Florida 32250







All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

UNICALL®

Retail Price List

Sept 1, 2011

Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
14001-06'	 <p>Classic Model Singles</p>	22.55	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14001-08'		24.70		
14001-10'		26.90		
14025-06'	 <p>Classic Model Duals</p>	32.30	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14025-08'		35.45		
14050-06'	 <p>Classic Model Singles</p>	36.30	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
14050-08'		38.15		
14050-10'		40.00		
14075-06'	 <p>Classic Model Duals</p>	48.54	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs std. and pillow speaker stations Edwards - Altec std. station Sylvania std. station Auth - All stations with no holding relays built in...	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
14075-08'		53.50		
14201-06'	 <p>Classic Model Singles</p>	22.65	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14201-08'		24.90		
14201-10'		27.35		
14225-06'	 <p>Classic Model Duals</p>	32.40	Executone, std-and pillow speaker stations. Dukane std (non-pillow speaker station) Picker Briggs - Edwards - Altec std. station Sylvania std. station and Auth std.	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14225-08'		35.80		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)









DWYER PRECISION® PRODUCTS INC.

266 N. 20th Street
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

UNICALL®		Retail Price List		Sept 1, 2011
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
14050-AU5-06'	 <p>Classic Model Singles</p>	64.01	Auth - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
14050-AU5-08'		66.24		
14050-AU5-10'		68.41		
14050-FA5-06'	 <p>Classic Model Singles</p>	64.01	Faraday - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
14050-FA5-08'		66.24		
14050-FA5-10'		68.42		
14050-ED5-06'	 <p>Classic Model Singles</p>	64.01	Edwards - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
14050-ED5-08'		66.24		
14050-ED5-10'		68.52		
14050-C05-06'	 <p>Classic Model Singles</p>	66.40	Couch - 5 - Pin Plug Maintain Contact	Normally open contacts...Pressing Bulb makes contact that is maintained until reset button (on top of plug at wall station) is pushed down. No system relay is required. Replaces pull chains or string actuated type. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC Suitable for use in oxygen atmosphere.
14050-C05-08'		68.57		
14050-C05-10'		70.79		
14007-ED3-06'	 <p>Classic Model Singles</p>	66.40	Edwards - 3 - Prong Plug Lifeline Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-ED3-08'		68.57		
14007-ED3-10'		70.80		
14007-WE8-06'	 <p>Classic Model Singles</p>	39.59	Intego / Wescom / Tektone / Simplex 8 pin connector	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-WE8-08'		41.76		
14007-WE8-10'		44.04		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)








DWYER PRECISION® PRODUCTS INC.

266 N. 20th Street
Jacksonville Beach, Florida 32250

All Nurse Call cords are suitable for use in oxygen atmosphere.

Phone#: (904) 249-3545 Fax#: (904) 249-1120

UNICALL®		Retail Price List		Sept 1, 2011
Catalog No.	Item Photo.	Price.	Systems Used In.	Electrical Characteristics
14007-EK8-06'	<div style="display: flex; justify-content: space-between;">Classic ModelSingles</div> 	66.40	Ektacom - 8 - Pin Octual In-Line Model with Keyway	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-EK8-08'		68.57		
14007-EK8-10'		70.79		
14007-EK15-06'	<div style="display: flex; justify-content: space-between;">Classic ModelSingles</div> 	55.79	Fisher Berkely - 15 - Pin Amphenol MMC / Lifeline In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-EK15-08'		58.01		
14007-EK15-10'		60.23		
14007-SY20-06'	<div style="display: flex; justify-content: space-between;">Classic ModelSingles</div> 	54.60	HCE- Sylvania / Phillips - 20 - Pin Plug In-Line Model	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-SY20-08'		56.87		
14007-SY20-10'		59.10		
14007-DN8-06'R	<div style="display: flex; justify-content: space-between;">Classic ModelSingles</div> 	40.05	8 - Pin Din Plug Rauland Borg System Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-DN8-08'R		42.28		
14007-DN8-10'R		44.60		
14007-DN8-06'J	<div style="display: flex; justify-content: space-between;">Classic ModelSingles</div> 	40.05	8 - Pin Din Plug Jeron Systems Standard Nurse Call	Momentary contact type. Normally open contacts, pressing bulb makes momentary contact and system relay holds call signal until reset button at bed or central nurse station is depressed. Electrical plug has tip and sleeve jack. Rated 1/3 amps 24 volt AC/DC
14007-DN8-08'J		42.28		
14007-DN8-10'J		44.60		

★ ★ ★ All prices are subject to change without notice, and are FOB DPPI ★ ★ ★

NOTE: (For classic celcon models other than what is listed add \$1.50 per unit).
Add \$1.00 per foot for longer length on cords. (Except Unicall or A & B Units - Not Available)